TiO MeshTM Surgical Mesh Implant

- DE Gebrauchsanweisung
- EN Instructions for Use
- FR Notice d'utilisation
- ES Instrucciones de uso
- PT Instruções de Uso
- IT Istruzioni per l'uso
- EL Οδηγίες χρήσης
- DA Brugsanvisning
- NL Gebruiksaanwijzing
- NO Bruksanvisning
- SV Bruksanvisning
- TR Kullanım Talimatları
- PL Instrukcja użytkowania
- CS Navod k použití
- SK Návod na použitie
- HU Használati utasítás
- RU Инструкция по применению

BioCer Entwicklungs-GmbH
 Ludwig-Thoma-Straße 36c
 95447 Bayreuth, Germany
 Tel.: +49 (0921) 787770-0
 Fax: +49 (0921) 787770-79
 info@biocer-gmbh.de
 www.biocer-gmbh.de

TiO MeshTM is not made with natural rubber latex.

Federal law (USA) restricts this device to sale by or on the order of physician.

FDA certified color additive ([Phthalocyaninato (2-)] Copper Blue) was used as a pigment in a concentration of max 0.5% by weight of the blue monofilament.

BioCer and TiO MeshTM are trademarks of Bio-Cer Entwicklungs-GmbH in Bayreuth.

This IFU is valid for both types, TiO MeshTM as well as TiO MeshTM light.

The link for the Summary of Safety and Clinical Performance (SSCP) can be found at https://www.biocer-gmbh.de/reports

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
ENGLISH

Description
 TiO MeshTM is made from a monofilament polypropylene thread and has a large-pore structure with blue orientation stripes. The surgical mesh implant has a titanium dioxide-coated surface. The titanium dioxide coating around the plastic has no influence on radiological diagnostic procedures (X-ray, CT, MRT and ultrasound scans). The fabric can be stretched in both directions and is highly flexible to react to body dynamics. TiO MeshTM is available in different shapes and sizes. Sterilization is performed with ethylene oxide.

Indications
 TiO MeshTM is a surgical mesh implant specially indicated for repair of tissue defects of the abdominal wall, where a non-reabsorbable support material is required. Relevant applications include the repair of inguinal and incisional hernia in all common surgical methods. TiO MeshTM can be used for both laparoscopic and other open procedures.

Contraindications
 Do not use TiO MeshTM in patients in the growth period, such as children and infants. The mesh implant is not for use in sites which are contaminated or infected.

Warnings
 The product is intended for single use only. Never re-use or re-sterilize any part of the product. Re-use of a single use product can lead to severe health problems and can cause the implant to malfunction. The implants are packaged in a sterile double-wrap packaging. Carefully check that the package is undamaged before use. If the package is damaged, open or moist, do not use the implant. Do not use damaged products. In order to guarantee the tensile strength of TiO MeshTM, contact with sharp objects should not exceed standard clinical practice during implantation in order to avoid damage to the mesh structure.

Side effects
 Possible adverse effects from the use of TiO MeshTM mesh implants include seroma, hematoma, recurrence, infection, visceral adhesions, mesh migration and fistula formation. Like all foreign bodies, this product may have a negative effect on pre-existing infection or contamination of the wound. As described in the literature, additional adverse effects including mesh shrinkage, chronic pain and reduction in sperm quality cannot be ruled out in connection with the use of mesh implants. Additional factors, such as age and gender, in addition to comorbidity and the extent of the hernia, may have an impact on adverse effects.


Precautions
 The product should only be used by doctors who are familiar with surgical techniques of this type. Before using TiO MeshTM, the physician should check the size of the hernial orifice and ensure a sufficient overlap. To minimize postoperative complications and recurrences, the mesh must be adequately affixed. This can be conducted in accordance with surgical practices using stitches (resorbable and non-resorbable), adhesives and staples. Take special caution when affixing the mesh with sutures, staples or tackers in the vicinity of nerves and vessels. When adjusting the meshes individually by cutting them, please bear in mind that by doing so you also remove the blunt laser-cut edges of the mesh which are intended to prevent microtraumata. In order to maintain the high tensile strength of TiO MeshTM, incisions and recesses in the direction of the blue strips must generally be avoided (does not apply to TiO MeshTM light). The manufacturer offers specially manufactured meshes for indications where an incision is necessary.

Post-operative recommendations
 Heavy lifting, physical work and sport should be avoided for 4-6 weeks.

Storage
 The recommended storage conditions are dry at room temperature, away from moisture and direct heat, protected against UV light, in the original packaging only, do not use past the expiration date.

Explanation of symbols on the packaging

Postoperative recommendations
 Heavy lifting, physical work and sport should be avoided for 4-6 weeks.

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