

HaemoCer™ - versatile use:

- Treatment of diffuse and localized bleeding
- For bleeding around suture lines
- Substitution of fibrin glue for haemostasis in vascular surgery
- Substitution of bone wax in bone bleeding
- Reduction in the need for electrocautery and associated tissue damage
- Use with compatible filter cell-savers

Fields of application:

- Cardiac, Vascular and Thoracic Surgery
- General Surgery
- Transplant Surgery
- Traumatology
- Tumour Surgery
- Oral and maxillofacial Surgery
- ENT Surgery
- Urology / Gynaecology
- Orthopaedics
- Plastic Surgery
- Burn surgery
- Neurosurgery
- Paediatrics

HaemoCer™

Haemostatic powder



Product overview HaemoCer™

HFP101	1g	5 pcs
HFP102	2g	5 pcs
HFP103	3g	5 pcs
HFP105	5g	5 pcs
HFZ101	440 mm	5 pcs
HFZ201	180 mm	5 pcs



100%
resorbable

made in
GERMANY

BioCer

BioCer Entwicklungs-GmbH
Ludwig-Thoma-Str. 36c · 95447 Bayreuth · Germany · Tel. +49 (0)921 78 77 70 0
Fax. +49 (0)921 78 77 70 79 · info@biocer-gmbh.de · www.biocer-gmbh.de



HDD103.10/2018-04

Natural - Safe - Effective

BioCer

HaemoCer™

HaemoCer™ Absorbable Polysaccharide Haemostat (APH) is a proprietary patent pending technology based on BioCers **Polysaccharide Ultra-hydrophilic Resorbable Engineering (PURE)** process. HaemoCer™ is a resorbable plant based haemostat containing no animal nor human components.

Haemostasis occurs by rapidly accelerating the normal physiologic clotting cascade without participation of chemical or pharmaceutical ingredients.

Natural

- Plant based
- Contains no animal or human components
- Complete resorption

Safe

- Rapid resorption minimizes risks of foreign body reactions and granuloma formation
- Reduction of homologous transfusions and blood donations

Effective

- Double enhanced absorption capacity, compact packaging reducing storage volume
- Rapid off-the-shelf deployment, no special storage or defrosting

Ease of use, DAPI application technique

Dry

Remove excess pooled blood by suction or gauze



Apply

Apply HaemoCer™ immediately and extensively to the wound site



Press

Apply gentle compression over the wound using dry gauze - in the case of profuse bleeding



Irrigate

Irrigate thoroughly prior to gently removing the gauze to ensure the gelled matrix remains at the wound site



Nature refined

HaemoCer™ – Absorbable Polysaccharide Haemostat (APH) engineering incorporates a sophisticated, plant-based polymer crosslinking that creates ultra-hydrophilic, biocompatible particles.

Once in contact with blood, HaemoCer™ enhances the natural clotting cascade by rapidly dehydrating the blood and accelerating the concentration of platelets, red blood cells and coagulation proteins at the bleeding site.

The second mode of action of HaemoCer™ is the formation of a robust gelled matrix that adheres to the bleeding site and forms a

mechanical barrier to prevent further bleeding. HaemoCer™ contains no animal nor human components and is rapidly resorbed by amylase within 48h without leaving any residues behind.

Latest studies demonstrate that plant based haemostatic materials (polysaccharides) like HaemoCer™ also reduce postoperative adhesions.

100%
resorbable

plant based

