

Endoscopic
Set

Made in
Germany

New
Device

HaemoCer™ Endoscopic SET

Haemostatic powder for
endoscopic application



Tube length 2400 mm



Outer tube diameter of 2.5 mm



HEPA filter -
High Efficiency Particulate Air




Usable instrument channel:
2.8 mm and above



BioCer

Endoscopic SET

Haemostatic Powder + Endoscopic Applicator



HaemoCer™ Haemostatic Powder

- Absorbable Polysaccharide Haemostat
- Containing no animal nor human components
- Plant-based polymer crosslinking that creates ultra-hydrophilic, biocompatible particles
- Rapidly resorbed by amylase within 48h without leaving any residues behind
- Rapid off-the-shelf deployment, no special storage, defrosting or premixing



HaemoCer™ Endoscopic Applicator

- For upper and lower gastrointestinal bleedings
- Even effective in difficult to access regions
- No blocking in case of contact with moisture
- Tube length 2400 mm
- Outer tube diameter of 2.5 mm
- HEPA filter - High Efficiency Particulate Air
- Usable instrument channel: 2.8 mm and above
- No high pressure gas source needed
- No compressor system needed
- Reduced dust formation
- Precise dosing of the applied powder

Product Overview HaemoCer™ Endoscopic Set

HFS323 HaemoCer™ PLUS 3 g (HFP203) + Endoscopic Applicator (HFZ103) 1 pcs

BioCer

BioCer Entwicklungs-GmbH

Ludwig-Thoma-Str. 36c · 95447 Bayreuth · Germany · Tel. +49 (0)921 78 77 70 0
Fax. +49 (0)921 78 77 70 79 · info@biocer-gmbh.de · www.biocer-gmbh.de



Management
System
EN ISO
13485:2016

www.tuv.com
ID 000059855

