

Universal
Applicator

made in
Germany

HaemoCerTM

Universal Applicator

For all minimally invasive surgeries

HaemoCerTM Universal Applicator is designed to deliver the haemostatic HaemoCerTM powder in laparoscopic, ENT, spinal and other MIS procedures. It guarantees precise application of APH particles to the bleeding site where this is mandatory. The delivery instrument can be easily attached to the HaemoCerTM or HaemoCerTM PLUS bellows dispenser and enables hemostat delivery under direction to the wound site for the control of capillary, venous and arteriole bleeding.

Fast - Effective - Safe



440 mm applicator



High flexibility



Integrated closure clip



Innovative applicator tip

BioCer



HaemoCer™

Universal Applicator



- Application in minimally invasive surgical procedures
- Longer applicator for deeper sourced bleeding
- Innovative applicator tip features a design which minimizes risks of blockage in case of contact with moisture
- One hand operated press-and-release locking system facilitating Universal Applicator tip placement without uncontrolled discharge
- Highly flexible enabling precise tip placement and offering enhanced handling and maximizing wound coverage capability
- Push-on coupling sheath offering simple and firm connection to HaemoCer™ or HaemoCer™ PLUS



Product Overview HaemoCer™ Universal Applicator

HFZ101

440 mm

5 pcs



BioCer

BioCer Entwicklungs-GmbH

Ludwig-Thoma-Str. 36c · 95447 Bayreuth · Germany · Tel. +49 (0)921 78 77 70 0
Fax. +49 (0)921 78 77 70 79 · info@biocer-gmbh.de · www.biocer-gmbh.de

