

FlexApplicator

made in
Germany

HaemoCer™ FlexApplicator

Fast - Effective - Safe

Malleable, efficient and varied

HaemoCer™ FlexApplicator is an accessory for our HaemoCer™ powder product line with standard applicator and can be installed easily. Because of its integrated wire, the HaemoCer™ FlexApplicator is bendable in different directions and still dimensionally stable.

HaemoCer™ or HaemoCer™ PLUS can be used at blind spots and complex challenging situations. HaemoCer™ FlexApplicator provides high flexibility and ensures a controlled and pinpoint dosage. A special rounded applicator tip prevents surgical lesions during application.



180 mm Applicator



Flexible





Dimensionally stable



Atraumatic
applicator tip

BioCer

HaemoCer™ FlexApplicator



- A simple plug mechanism for fast and firm connection to our HaemoCer™ powder product line with standard applicator
- Longer applicator for bleeding sites that are difficult to access
- HaemoCer™ FlexApplicator is bendable at several points at the same time and dimensionally stable to reach even difficult surgical sites.
- The transparency of the applicator enables a clear view to the amount of HaemoCer™ powder applied, therefore a controlled and pinpointed application is guaranteed
- Even after multiple bendings the powder freely flows through the applicator
- Atraumatic applicator tip prevents surgical lesions

Product Overview HaemoCer™ FlexApplicator

HFZ201

180 mm

5 pcs

BioCer

BioCer Entwicklungs-GmbH

Ludwig-Thoma-Str. 36c · 95447 Bayreuth · Germany · Tel. +49 (0)921 78 77 70 0

Fax. +49 (0)921 78 77 70 79 · info@biocer-gmbh.de · www.biocer-gmbh.de

