

Flex Set

made in
Germany

HaemoCer™ Flex SET

Fast - Effective - Safe



180 mm Applicator



Flexible



Dimensionally stable




Atraumatic
applicator tip



BioCer

Flex SET

Haemostatic Powder + FlexApplicator
Malleable, efficient and varied



Haemostatic Powder

- Treatment of diffuse and localized bleeding
- For bleeding around suture lines
- Substitution of fibrin glue for haemostasis in vascular surgery
- Substitution of bone wax in bone bleeding
- Reduction in the need for electrocautery and associated tissue damage
- Use with compatible filter cell-savers



HaemoCer™ FlexApplicator

- Longer applicator for bleeding sites that are difficult to access
- HaemoCer™ FlexApplicator is bendable at several points at the same time and dimensionally stable to reach even difficult surgical sites
- The transparency of the applicator enables a clear view to the amount of HaemoCer™ powder applied, therefore a controlled and pinpointed application is guaranteed
- Even after multiple bendings the powder freely flows through the applicator
- Atraumatic applicator tip prevents surgical lesions

Product Overview HaemoCer™ Flex Set

HFS215	HaemoCer™ 5g (HFP105) + 180 mm Applicator (HFZ201)	1 pcs
HFS225	HaemoCer™ PLUS 5g (HFP205) + 180 mm Applicator (HFZ201)	1 pcs

BioCer

BioCer Entwicklungs-GmbH

Ludwig-Thoma-Str. 36c · 95447 Bayreuth · Germany · Tel. +49 (0)921 78 77 70 0
Fax. +49 (0)921 78 77 70 79 · info@biocer-gmbh.de · www.biocer-gmbh.de

