

Universal
Set

made in
Germany

HaemoCer™ Universal SET

Fast - Effective - Safe



440 mm Applicator



High flexibility



Integrated closure clip




Innovative applicator tip



BioCer

Universal SET

Haemostatic Powder + Universal Applicator
For all minimally invasive surgeries



Haemostatic Powder

- Treatment of diffuse and localized bleeding
- For bleeding around suture lines
- Substitution of fibrin glue for haemostasis in vascular surgery
- Substitution of bone wax in bone bleeding
- Reduction in the need for electrocautery and associated tissue damage
- Use with compatible filter cell-savers



HaemoCer™ Universal Applicator

- Longer applicator for deeper sourced bleeding
- Innovative applicator tip features a design which minimizes risks of blockage in case of contact with moisture
- One hand operated press-and-release locking system facilitating Universal Applicator tip placement without uncontrolled discharge
- Highly flexible enabling precise tip placement and offering enhanced handling and maximizing wound coverage capability
- Push-on coupling sheath offering simple and firm connection to HaemoCer™ or HaemoCer™ PLUS

Product Overview HaemoCer™ Universal Set

HFS115	HaemoCer™ 5g (HFP105) + 440 mm Applicator (HFZ101)	1 pcs
HFS125	HaemoCer™ PLUS 5g (HFP205) + 440 mm Applicator (HFZ101)	1 pcs

BioCer

BioCer Entwicklungs-GmbH

Ludwig-Thoma-Str. 36c · 95447 Bayreuth · Germany · Tel. +49 (0)921 78 77 70 0
Fax. +49 (0)921 78 77 70 79 · info@biocer-gmbh.de · www.biocer-gmbh.de

